

GUIDE TO *Product* USE



GUIDE TO PRODUCT USE

Table of Contents

1	Anaerobes	8	<i>Legionella</i> spp.
1	<i>Acanthamoeba</i>	8	<i>Listeria</i> spp.
1	<i>Acinetobacter baumannii</i>	9	<i>Micrococcus</i> spp.
2	<i>Aeromonas hydrophila</i>	9	<i>Moraxella catarrhalis</i>
2	<i>Bacillus</i> spp.	9	<i>Mycobacterium</i> spp.
2	<i>Bordetella pertussis</i>	9	<i>Mycoplasma</i> spp.
2	<i>Borrelia burgdorferi</i> (lyme disease)	10	<i>Naegleria</i> spp.
2	<i>Brucella</i> spp.	10	<i>Neisseria gonorrhoeae</i>
2	<i>Burkholderia cepacia</i>	10	<i>Nocardia</i> spp.
2	<i>Campylobacter jejuni</i>	10	Parasitology
3	<i>Candida</i> spp.	10	<i>Plasmodium</i> spp. (malaria)
3	<i>Chlamydia</i> spp.	10	<i>Pneumocystis</i>
3	<i>Chronobacter sakazakii</i>	10	<i>Pseudomonas aeruginosa</i>
3	<i>Clostridium difficile</i>	11	<i>Salmonella</i> spp.
3	<i>Coccidioides immitis</i>	11	<i>Shigella</i> spp.
3	Coliforms	12	<i>Staphylococcus aureus</i>
4	<i>Corynebacterium</i> spp.	12	<i>Staphylococcus aureus</i> , Methicillin Resistant
4	<i>Cryptococcus</i> spp.	12	<i>Streptococcus</i> group A
4	<i>Cryptosporidium</i> spp.	12	<i>Streptococcus</i> group B
4	Dermatophytic Fungi	12	<i>Streptococcus anginosus</i>
5	<i>Enterobacteriaceae</i> (also see Coliforms)	12	<i>Streptococcus</i> (Beta-hemolytic) A, B, C, D, F, G
6	<i>E. coli</i> and <i>E. coli</i> O157 (also see Coliforms)	13	<i>Streptococcus pneumoniae</i>
6	<i>Enterococcus</i> spp. and VRE	13	<i>Toxoplasmosis</i> spp.
7	Environmental Monitoring of Surfaces	13	<i>Treponema pallidum</i> (syphilis)
7	Fungi, Molds	13	<i>Trichomonas</i> spp.
7	<i>Gardnerella vaginalis</i>	13	<i>Ureaplasma</i>
7	<i>Giardia lamblia</i>	13	<i>Vibrio</i> spp.
8	<i>Haemophilus</i> spp.	13	<i>Yersinia enterocolitica</i>
8	<i>Helicobacter pylori</i>	13	Virus
8	<i>Lactobacillus</i> spp.	13	Miscellaneous



Headquarters

1430 West McCoy Lane
Santa Maria, CA 93455
800.266.2222 : phone
805.346.2760 : fax
Sales@HardyDiagnostics.com
www.HardyDiagnostics.com

Distribution Centers

Santa Maria, California
Centralia, Washington
Salt Lake City, Utah
Phoenix, Arizona
Dallas, Texas
Springboro, Ohio
Lake City, Florida
Albany, New York



ISO 13485
FM 572526

The Quality Management System at the Hardy Diagnostics manufacturing facility is certified to ISO 13485.

Copyright © 2011 Hardy Diagnostics

Guide to Product Use



Organism	Usage	Product	Cat. no.	Section	Page
Anaerobes	Cultivation and Isolation	<i>See Anaerobic Microbiology for pre-reduced media</i>			
		AnaBag 150™	AB20	Anaerobic Microbiology	A7
		Anaerobe Jar, Complete	16000	Anaerobic Microbiology	A12
		Anaerobic Indicator	BR55	Anaerobic Microbiology	A6
		Anaerobic PEA	A90	Prepared Culture Media	N1
		AnaeroGen™ Anaerobic Gas Generator	AN25US	Anaerobic Microbiology	A7
		AnaeroGen™ Compact	AN010C	Anaerobic Microbiology	A6
		AnaeroGRO™ Products	MULTIPLE	Anaerobic Microbiology	A1-A4
		BBE Agar	G05	Prepared Culture Media	N3
		BBE/LKV Biplate	J102	Prepared Culture Media	N3
		Brucella Agar with Hemin and Vitamin K/LKV	J87	Prepared Culture Media	N8
		Brucella Blood Agar with Hemin and Vitamin K	A30	Prepared Culture Media	N8
		Brucella with H&K/Anaerobic PEA	J74	Prepared Culture Media	N8
		LKV Agar	A60	Prepared Culture Media	N26
		Neomycin Anaerobic Blood Agar	A62	Prepared Culture Media	N33
		PEA/LKV Biplate	J86	Prepared Culture Media	N35
	Thioglycollate with H/K and Fluid Thioglycollate	MULTIPLE	Prepared Culture Media	—	
	Identification	Anaerobe Nitrate Reagent A	Z134	Anaerobic Microbiology	A13
		Anaerobe Nitrate Reagent B	Z135	Anaerobic Microbiology	A13
		Bile HardyDisks™	Z7091	Anaerobic Microbiology	A13
		Catalase Reagent	Z262	Anaerobic Microbiology	A13
		Colistin Sulfate, CT, HardyDisk™	Z8411	Anaerobic Microbiology	A14
		Egg Yolk Agar, Modified	G215	Prepared Culture Media	N16
		Esculin Agar Slant	L29	Prepared Culture Media	N17
		Indole Nitrate Broth	K147	Prepared Culture Media	N23
		Kanamycin	Z7191	Rapid Tests & Reagents	P5
		Nitrate, HardyDisk™	Z7071	Rapid Tests & Reagents	P8
		Nitrate Reagent C	Z73	Rapid Tests & Reagents	P1
		Nitrocef Matchbook™	Z108	Rapid Tests & Reagents	P9
		PRO Disks	K1532	Anaerobic Microbiology	A15
		Rapid Urea/Esculin Kit	Z57	Rapid Tests & Reagents	P11
		SPS Disks	Z7381	Rapid Tests & Reagents	P12
	Vancomycin	Z7501	Rapid Tests & Reagents	P14	
	Transport	Anaerobic Media, without Swab	S1200	Anaerobic Microbiology	A5
		Anaerobic Media, with Rayon Swabs	S120D	Anaerobic Microbiology	A5
		Anaerobic Transport Medium	AG25H	Anaerobic Microbiology	A5
	Quality Control	<i>B. fragilis</i> , ATCC® 25285™, KWIK-STIK™ DuoPak	0320P	Control Organisms	D4
		<i>B. thetaiotaomicron</i> , ATCC® 29741™, KWIK-STIK™ DuoPak	0319P	Control Organisms	D4
		<i>C. difficile</i> , ATCC® 9689™, KWIK-STIK™ DuoPak	0329P	Control Organisms	D6
		<i>C. perfringens</i> , ATCC® 13124™, KWIK-STIK™ DuoPak	0318P	Control Organisms	D6
		<i>E. lenta</i> , ATCC® 43055™, KWIK-STIK™ DuoPak	0936P	Control Organisms	D7
Acanthamoeba	Cultivation and Isolation	Non-Nutrient Agar	G225	Prepared Culture Media	N33
		Page's Saline	R225	Prepared Culture Media	N34
	Stain	Modified Field's Stain Kit	Z15	Stains & Reagents	Q16
Maintenance of Stock Cultures	Acanthamoeba Broth (PYG)	K225	Prepared Culture Media	N1	
Acinetobacter baumannii	Cultivation and Isolation	Acinetobacter Biplate	J135	Prepared Culture Media	N1
		Leeds Medium	G261	Prepared Culture Media	N25
		MDR Acinetobacter Medium	G259	Prepared Culture Media	N29
	Quality Control	<i>A. baumannii</i> , ATCC® BAA-1605™, KWIK-STIK™ DuoPak	01057P	Control Organisms	D2
		<i>A. baumannii</i> , ATCC® BAA-747™, KWIK-STIK™ DuoPak	0119P	Control Organisms	D2
		<i>B. cepaciae</i> , ATCC® 25416™, KWIK-STIK™ DuoPak	0488P	Control Organisms	D4
		<i>E. coli</i> , ATCC® 25922™, KWIK-STIK™ DuoPak	0335P	Control Organisms	D8

Guide to Product Use



Organism	Usage	Product	Cat. no.	Section	Page
<i>Aeromonas hydrophila</i>	Cultivation and Isolation	Ampicillin Dextrin Agar	G180	Prepared Culture Media	N1
		Blood Agar with Ampicillin	A12	Prepared Culture Media	N6
		CIN Agar/MacConkey Agar with Sorbitol	J49	Prepared Culture Media	N12
		CIN Agar	G20	Prepared Culture Media	N12
Quality Control	A. hydrophila, ATCC® 7966™, KWIK-STIK™ DuoPak	0870P	Control Organisms	D2	
		E. coli, ATCC® 25922™, KWIK-STIK™ DuoPak	0335P	Control Organisms	D8
<i>Bacillus</i> spp.	Cultivation and Isolation	AK #2 Agar	G67, U189	Prepared Culture Media	N1
		Bacillus cereus Medium, CRITERION™	C5061	Dehydrated Culture Media	G6
		Bacillus cereus Selective Agar	G156	Prepared Culture Media	N2
		MYP Agar	G147	Prepared Culture Media	N32
		Tryptic Soy Broth (TSB) with Polymyxin B	K180	Prepared Culture Media	N50
	Identification	CRET-RPLA, B. cereus Enterotoxin Detection Kit Microgen® Bacillus ID Panel	230102	Direct Specimen Testing	H2
MID66			Rapid Tests & Reagents	P7	
<i>Bordetella pertussis</i>	Cultivation and Isolation	Bordet Gengou Agar	A05	Prepared Culture Media	N7
		Regan-Lowe Agar	A65	Prepared Culture Media	N39
		Regan-Lowe Agar without Cephalexin	A63	Prepared Culture Media	N39
		Regan-Lowe Semi Solid Agar Deep	Q32	Prepared Culture Media	N39
	Transport	Regan-Lowe Semi Solid Agar Deep	Q32	Prepared Culture Media	N39
	Quality Control	B. pertussis, ATCC® 9340™, KWIK-STIK™ DuoPak B. pertussis, ATCC® 9797™, KWIK-STIK™ DuoPak	0100P	Control Organisms	D4
0843P			Control Organisms	D4	
<i>Borrelia burgdorferi</i> (lyme disease)	Identification	ImmunoDot™, Borrelia	2470445	Direct Specimen Testing	H10
		ImmunoDot™, Borrelia, EIA	2470882	Direct Specimen Testing	H10
		ImmunoWell®, Borrelia, EIA	2470437	Direct Specimen Testing	H10
		Premier™, Human Lyme	696016	Direct Specimen Testing	H10
<i>Brucella</i> spp.	Identification	Basic Fuchsin, 1:25,000 dilution	Z005	Prepared Culture Media	N8
		Basic Fuchsin, 1:50,000 dilution	Z006	Prepared Culture Media	N8
		Basic Fuchsin, 1:100,000 dilution	Z007	Prepared Culture Media	N8
		Thionin Acetate, 1:25,000 dilution	Z035	Prepared Culture Media	N8
		Thionin Acetate, 1:50,000 dilution	Z036	Prepared Culture Media	N8
		Thionin Acetate, 1:100,000 dilution	Z037	Prepared Culture Media	N8
	Serology	Febrile Antigen Set	240786	Antisera	B14
<i>Burkholderia cepacia</i>	Cultivation and Isolation	BCSA	G09	Prepared Culture Media	N3
		PC Agar	G48	Prepared Culture Media	N35
	Quality Control	B. cepaciae, ATCC® 25416™, KWIK-STIK™ DuoPak	0488P	Control Organisms	D4
<i>Campylobacter jejuni</i>	Enrichment	Campy Thioglycollate	K128	Prepared Culture Media	N10
	Cultivation and Isolation	AnaeroGRO™ Campylobacter Selective Agar	AG701	Anaerobic Microbiology	A2
		Bolton Broth	U83	Prepared Culture Media	N7
		Campy, Blood Free, Karmali Agar	G06	Prepared Culture Media	N10
		Campy Cefex Agar, Modified	A122	Prepared Culture Media	N10
		Campy CVA Agar with 5% Sheep Blood	A40	Prepared Culture Media	N10
		CampyGen Gas Generator	CN025A	Anaerobic Microbiology	A9
	Identification	Campylobacter Latex Test, Microgen® Hippurate Test Indoxyl Acetate Disks Ninhydrin Reagent	F46	Rapid Tests & Reagents	P2
			Z52	Rapid Tests & Reagents	P4
			Z111	Rapid Tests & Reagents	P5
			Z70	Rapid Tests & Reagents	P8
Quality Control	C. jejuni, ATCC® 33291™, KWIK-STIK™ DuoPak	0481P	Control Organisms	D4	

Guide to Product Use



Organism	Usage	Product	Cat. no.	Section	Page
Candida spp.	Cultivation and Isolation	BiGGY Agar	G17, L14	Prepared Culture Media	N5
		Candida Isolation Agar	G16	Prepared Culture Media	N11
		Emmon's Sabdex Agar	W20	Prepared Culture Media	N40
		HardyCHROM™ Candida	G301	Prepared Culture Media	N20
		Mycobiotic Agar	L45, W50, X30	Prepared Culture Media	N32
		Sabdex Agar	W70	Prepared Culture Media	N40
	Identification	AlbiQuick™	Z121	Rapid Tests & Reagents	P1
		Corn Meal Agar with Tween® 80	L55, Q12, W10	Prepared Culture Media	N13
		Germ Tube Cryo™	Z217	Rapid Tests & Reagents	P4
		GlabrataQuick™	Z298	Rapid Tests & Reagents	P4
		Rapid Trehalose Broth	Z205	Prepared Culture Media	N47
		Trehalose Fermentation Broth	Y205	Prepared Culture Media	N47
Quality Control	<i>C. glabrata</i> , ATCC® 2001™, KWIK-STIK™ DuoPak	0992P	Control Organisms	D5	
	<i>C. albicans</i> , ATCC® 10231™, KWIK-STIK™ DuoPak	0443P	Control Organisms	D5	
	<i>C. glabrata</i> , ATCC® 66032™, KWIK-STIK™ DuoPak	0986P	Control Organisms	D5	
Chlamydia spp.	Diagnostic Test	QuickVue® Chlamydia	B006	Direct Specimen Testing	H2
	Transport	Clearview™ Chlamydia Collection Kit	135315	Direct Specimen Testing	H2
		CVM Flocked Swab Kit Female	R96KF	Direct Specimen Testing	H2
		CVM Flocked Swab Kit Male	R96KM	Direct Specimen Testing	H2
		CVM Transport Media	R96	Direct Specimen Testing	H2
Chronobacter sakazakii	Cultivation and Isolation	HardyCHROM™ Sakazakii	G315	Prepared Culture Media	N21
	Quality Control	<i>C. mytjensii</i> , ATCC® 51329™, KWIK-STIK™ DuoPak	0756P	Control Organisms	D7
Clostridium difficile	Cultivation and Isolation	AnaeroGRO™ CCFA	AG501	Anaerobic Microbiology	A2
	Diagnostic Test	<i>C. Diff</i> Quik Chek Complete®	30525C	Direct Specimen Testing	H3
		Clearview™ Clostridium difficile Toxin A	135050	Direct Specimen Testing	H3
		ImmunoCard® <i>C. difficile</i>	706050	Direct Specimen Testing	H3
		Premier™ Toxins A and B Enzyme Immunoassay	616096	Direct Specimen Testing	H4
Identification	PRO Disks	K1532	Anaerobic Microbiology	A15	
Quality Control	<i>C. difficile</i> , ATCC® 9689™, KWIK-STIK™ DuoPak	0329P	Control Organisms	D6	
Coccidioides immitis	Cultivation and Isolation	BHI Agar with Blood, Chloramphenicol and Gentamicin Slant	L114	Prepared Culture Media	N4
		BHI Agar with Blood and Gentamicin	X12	Prepared Culture Media	N4
		BHI Agar with Blood Slant	L31	Prepared Culture Media	N4
	Diagnostic Test	Premier™, Coccidioides	603096	Direct Specimen Testing	H4
Coliforms	Cultivation and Isolation	HardyCHROM™ ECC	G303	Prepared Culture Media	N21
		Violet Red Bile Agar with Glucose	G178	Prepared Culture Media	N52
	Identification	Microgen® Gram Negative ID, Panel A	MID64	Rapid Tests & Reagents	P7
	Water	Brilliant Green Bile Broth	K09, K66	Prepared Culture Media	N7
		Colitag™ P/A Test	46000013	Featured Products	T10
		EC Broth	K13, K14	Prepared Culture Media	N16
		EC Broth with MUG	K18	Prepared Culture Media	N16
		Lauryl Sulfate with MUG	K35	Prepared Culture Media	N25
		Lauryl Tryptose Broth	K32, K33	Prepared Culture Media	N25
		m Endo LES Agar	G28, G128, Q17	Prepared Culture Media	N27
m FC Media with 1% Rosolic Acid		G126, G127	Prepared Culture Media	N27	
Quality Control	<i>E. coli</i> , ATCC® 25922™, KWIK-STIK™ DuoPak	0335P	Control Organisms	D8	
	<i>K. pneumoniae</i> , ATCC® 31488™, KWIK-STIK™ DuoPak	0261P	Control Organisms	D10	
	<i>P. aeruginosa</i> , ATCC® 10145™, KWIK-STIK™ DuoPak	0416P	Control Organisms	D14	
	<i>E. aerogenes</i> , ATCC® 13048™, KWIK-STIK™ DuoPak	0306P	Control Organisms	D19	

Guide to Product Use




Organism	Usage	Product	Cat. no.	Section	Page	
Corynebacterium spp.	Cultivation and Isolation	Aerobic Corynebacterium Phosphomycin Agar	G162	Prepared Culture Media	N1	
		Loeffler's Slant	L28	Prepared Culture Media	N27	
		Serum Tellurite Agar Slant	L98	Prepared Culture Media	N42	
		Tinsdale Base	Q56	Prepared Culture Media	N46	
	Identification	RSU Broth	Z53	Rapid Tests & Reagents	P12	
Staining	Methylene Blue Stain	Z88	Stains & Reagents	Q14		
Quality Control	<i>C. diphtheriae</i> , ATCC® 13812™, KWIK-STIK™ DuoPak	0844P	Control Organisms	D6		
Cryptococcus spp.	Cultivation and Isolation	Caffeic Acid Agar Slant	L13	Prepared Culture Media	N10	
		CGB Agar	G113	Prepared Culture Media	N11	
	Diagnostic Test	Cryptococcus Antigen Latex Agglutination Premier™, Cryptococcal Antigen	140100 602096	Direct Specimen Testing Direct Specimen Testing	N5 N5	
		Identification	Caffeic Acid Disks	Z118	Rapid Tests & Reagents	P2
	Rapid Urea Medium		Z54	Rapid Tests & Reagents	P11	
	Urea Disks		DU7525	Rapid Tests & Reagents	P14	
	Staining	India Ink	Z64	Stains & Reagents	Q9	
Quality Control	<i>C. neoformans</i> , ATCC® 14116™, KWIK-STIK™ DuoPak	0291P	Control Organisms	D7		
	<i>C. gattii</i> , ATCC® MYA-4560™, KWIK-STIK™ DuoPak	01051P	Control Organisms	D7		
Cryptosporidium spp.	Diagnostic Test	ImmunoCard STAT!®, Cryptosporidium/Giardia	750830	Direct Specimen Testing	H6	
		Merifluor® Crypto/Giardia	250050	Direct Specimen Testing	H6	
		Para-Tect™ Cryptosporidium Antigen Detection Kit	MCCCP96	Direct Specimen Testing	H6	
		Para-Tect™ Cryptosporidium/Giardia Direct Fluorescent Assay	MCCCG75	Direct Specimen Testing	H6	
	Staining	Carbol Fuchsin for AFB	CF008	Stains & Reagents	Q1	
		Cryptosporidium Stain, Modified (Methylene Blue)	482K	Parasitology	M7	
Cryptosporidium Stain, Modified (Methylene Green)	482KMG	Parasitology	M7			
Quality Control	Control Slide, Cryptosporidium	SL4510	Parasitology	M10		
Dermatophytic Fungi	Cultivation and Isolation	Derm-Duet™ II, DTM/ESA	DMT10	Prepared Culture Media	N15	
		Derm-Duet™ RSM/DTM	J350	Prepared Culture Media	N15	
		DTM Agar Slant	L27, L115, X15	Prepared Culture Media	N16	
		Mycobiotic Agar	L45, W50, X30	Prepared Culture Media	N32	
		Mycovue™ Slide Culture System, Corn Meal Agar with Tween® 80	MV3	Prepared Culture Media	N32	
		Mycovue™ Slide Culture System, Potato Flake Agar	MV1	Prepared Culture Media	N32	
		Mycovue™ Slide Culture System, SABHI® Agar with Chloramphenicol and Gentamicin	MV2	Prepared Culture Media	N32	
		Potato Flake Selective Agar with Ciprofloxacin and Vancomycin	W159	Prepared Culture Media	N38	
		Sabdex Agar with Chloramphenicol	W72	Prepared Culture Media	N41	
		Identification	Dermatophyte Milk Agar	L75	Prepared Culture Media	N15
			Potato Dextrose Agar	L90, Q29	Prepared Culture Media	N37
			Potato Flake Agar	W59	Prepared Culture Media	N37
			Trichophyton #1	C51	Prepared Culture Media	N47
	Trichophyton #2		C52	Prepared Culture Media	N47	
	Trichophyton #3		C53	Prepared Culture Media	N47	
	Trichophyton #4		C54	Prepared Culture Media	N47	
	Trichophyton #5		C55	Prepared Culture Media	N47	
	Trichophyton #6	C56	Prepared Culture Media	N47		
	Trichophyton #7	C57	Prepared Culture Media	N47		
	Staining, Wet Mounts	Lactophenol Cotton Blue	Z68	Stains & Reagents	Q9	
		SoluPhyte™	Z178	Rapid Tests & Reagents	P12	

Guide to Product Use



Organism	Usage	Product	Cat. no.	Section	Page
Enterobacteriaceae (also see Coliforms)	Enrichment	EE Broth Mossel	K191, U291, U391	Prepared Culture Media	N16
	Cultivation and Isolation	CIN Agar for <i>Yersinia</i>	G20	Prepared Culture Media	N12
		CIN Agar/MacConkey Agar with Sorbitol	J49	Prepared Culture Media	N12
		EMB Levine Agar	G25	Prepared Culture Media	N16
		HardyCHROM™ Carbapenemase	G323	Prepared Culture Media	N20
		HardyCHROM™ ESBL	G321	Prepared Culture Media	N21
		HardyCHROM™ Urine Biplate	J100	Prepared Culture Media	N21
		HardyCHROM™ UTI	G313	Prepared Culture Media	N22
		MacConkey Agar	G35	Prepared Culture Media	N28
		TCBS Agar for <i>Vibrio</i>	G55, Q52, R23	Prepared Culture Media	N45
	Identification	Enteroscreen 4™	L225	Prepared Culture Media	N16
		Indole, DMACA Spot	Z65	Rapid Tests & Reagents	P5
		Koser Citrate Broth	R135	Prepared Culture Media	N24
		Microgen® GN- ID Panel A	MID66	Rapid Tests & Reagents	P7
		Microgen® GN- ID Panel B	MID65	Rapid Tests & Reagents	P7
		MIL Medium	Q09	Prepared Culture Media	N30
		MIO Medium	R24	Prepared Culture Media	N30
		MOT Disks	K9101	Rapid Tests & Reagents	P8
		MUG Disks	DM9525	Rapid Tests & Reagents	P8
		ONPG Broth	L05	Prepared Culture Media	N34
		ONPG Disks	Z7321	Rapid Tests & Reagents	P9
		ONPG-PDA Broth	R93	Prepared Culture Media	N34
		Phenylalanine Agar Slant	L21	Prepared Culture Media	N36
		Rapid Urea Medium	Z54, R154	Prepared Culture Media	N51
		Simmons Citrate Agar Slant	L80	Prepared Culture Media	N43
		TSI Agar Slant	R32	Prepared Culture Media	N47
		Urea Agar Slant	L65, R42	Prepared Culture Media	N51
		Urea Disks	DU7525	Rapid Tests & Reagents	P14
	Urease Test Tablets	K650	Rapid Tests & Reagents	P14	
	Quality Control	<i>E. coli</i> , ATCC® 25922™, KWIK-STIK™ DuoPak	0335P	Control Organisms	D8
		<i>E. coli</i> , ATCC® 35218™, KWIK-STIK™ DuoPak	0495P	Control Organisms	D8
		<i>K. pneumoniae</i> , ATCC® 700603™, KWIK-STIK™ DuoPak	0784P	Control Organisms	D10
		<i>P. mirabilis</i> , ATCC® 12453™, KWIK-STIK™ DuoPak	0440P	Control Organisms	D14



...an encyclopedia.

HUGO™ contains over 4,000 pages of reference material for the microbiologist:

- A complete description of each genus of medically important bacteria and fungi.
- A complete glossary of microbiology terms; the most extensive ever assembled.
- Includes parasitology, mycology, bacteriology, and epidemiology terms.
- Abbreviations included.
- Current taxonomy changes.
- A directory of organizations and agencies that relate to the microbiology laboratory, including proficiency testing programs. Phone numbers and websites included.
- Drug of choice listed with each pathogen.
- Organisms listed by body site with the degree of pathogenicity.
- Information on stains, dyes, and pH indicators used in microbiology.
- Reportable infectious diseases and their incidence in the U.S.
- Normal reference values of laboratory tests.
- Explanations of statistical terms used in QC.
- List of common stock organisms and how to efficiently maintain them.
- Information on microbiology stains and control slides for QC.
- Suggested antimicrobial disks to use when testing clinical isolates.
- Chart of recommended media to use when setting up incoming specimens.
- Answers to frequently asked questions regarding culture media.

Guide to Product Use



Organism	Usage	Product	Cat. no.	Section	Page
<i>E. coli</i> and <i>E. coli</i> O157 (also see Coliforms)	Enrichment	MacConkey Broth	K194, U125	Prepared Culture Media	N28
	Cultivation and Isolation	CNA/BluEcoli™ Urine Biplate	J116	Prepared Culture Media	N20
		BluEcoli™ Urine Biplate	J123	Prepared Culture Media	N20
		CT-SMAC with Cefixime and Tellurite (O157)	G129	Prepared Culture Media	N14
		EMB Levine Agar	G25	Prepared Culture Media	N16
		HardyCHROM™ ECC	G303	Prepared Culture Media	N21
		HardyCHROM™ <i>E. coli</i> O157 (food only)	G305	Prepared Culture Media	N21
		HardyCHROM™ Urine Biplate	J100	Prepared Culture Media	N21
		MacConkey Agar	G35	Prepared Culture Media	N28
		MacConkey Agar with Sorbitol (SMAC) (O157)	G36	Prepared Culture Media	N28
	Nutrient Agar with MUG	G114	Prepared Culture Media	N34	
	Identification	ColiScreen™	Z104	Rapid Tests & Reagents	P3
		<i>E. coli</i> PRO™ O157 Latex Kit	PL070HD	Rapid Tests & Reagents	P3
		Indole, DMACA Spot	Z65	Rapid Tests & Reagents	P5
		MacConkey Broth	U125	Prepared Culture Media	N28
		MOT disks	K9101	Rapid Tests & Reagents	P8
		MUG Disks	DM9525	Rapid Tests & Reagents	P8
		ONPG Broth	L05	Prepared Culture Media	N34
		ONPG Disks	Z7321	Rapid Tests & Reagents	P9
	Water	Brilliant Green Bile Broth	K09, K66	Prepared Culture Media	N7
		Colitag™ P/A Test	46000013	Featured Products	T10
		EC Broth, Double Strength	K14	Prepared Culture Media	N16
		EC Broth with MUG	K18	Prepared Culture Media	N16
		EC Broth	K13	Prepared Culture Media	N16
		Lauryl Sulfate with MUG	K35	Prepared Culture Media	N25
		Lauryl Tryptose Broth	K32, K33	Prepared Culture Media	N25
	Quality Control	<i>E. coli</i> , ATCC® 25922™, KWIK-STIK™ DuoPak	0335P	Control Organisms	D19
		<i>E. coli</i> (O157:H7), ATCC® 35150™, KWIK-STIK™ DuoPak	0617P	Control Organisms	D8
<i>E. coli</i> ATCC® 8739™, KWIK-STIK™ DuoPak		0483P	Control Organisms	D8	
<i>Enterococcus</i> spp. and VRE	Cultivation and Isolation	BE with Azide Broth	K11	Prepared Culture Media	N4
		BE Agar	G12, L10	Prepared Culture Media	N3
		BE Agar with Azide	G11, J64	Prepared Culture Media	N4
		BE Agar with Azide /CNA	J66	Prepared Culture Media	N4
		BE Agar with Azide, Vancomycin, Gentamicin, and Amphotericin B	G107	Prepared Culture Media	N4
		BE with Azide Broth with Vancomycin	R98	Prepared Culture Media	N4
		BE with Azide and Vancomycin	G103	Prepared Culture Media	N3
		Blood Agar, 5%	A10	Prepared Culture Media	N6
		Columbia CNA Agar	A50	Prepared Culture Media	N12
		Blood Agar, 5%/Bile Esculin Agar Biplate	J38	Prepared Culture Media	N5
	HardyCHROM™ VRE	G333	Prepared Culture Media	N22	
	Identification	LAP Kit	DL10K	Rapid Tests & Reagents	P5
		PYR Test Kit	Z75	Rapid Tests & Reagents	P10
		Rapid MGP Media	Z225	Rapid Tests & Reagents	P7
		StrepQuick™	Z122	Rapid Tests & Reagents	P13
	Susceptibility	BHI Agar with Vancomycin	G14	Prepared Culture Media	N5
	Environment	HardyCHROM™ VRE Agar Contact Plate	P19	Prepared Culture Media	N23
	Water	Azide Dextrose Broth	K86, K87	Prepared Culture Media	N2
		m EI Agar	G124	Prepared Culture Media	N27
	Quality Control	<i>E. faecalis</i> , ATCC® 29212™, KWIK-STIK™ DuoPak	0366P	Control Organisms	D22
		<i>E. faecalis</i> , ATCC® 51299™, KWIK-STIK™ DuoPak	0959P	Control Organisms	D8

Guide to Product Use



Organism	Usage	Product	Cat. no.	Section	Page	
Environmental Monitoring of Surfaces	Cultivation and Isolation	Baird Parker Agar	P16	Prepared Culture Media	N2	
		Blood Agar, 5%	P33	Prepared Culture Media	N6	
		Cellulose Agar	P83	Prepared Culture Media	N11	
		D/E Neutralizing Agar	P99	Prepared Culture Media	N14	
		D/E Neutralizing Agar, Irradiated, Triple Bagged	P13	Prepared Culture Media	N14	
		EnviroBootie™, EnviroBootie™ II	EB20, EB20M, EB60K	Swabs & Transports	S15	
		EnviroTrans™ Swab Rinse Kit Lethen Broth	SRK25	Prepared Culture Media	N17	
		EnviroTrans™ Swab Rinse Kit Neutralizing Buffer	SRK15	Prepared Culture Media	N17	
		HardyCHROM™ ECC	P17	Prepared Culture Media	N21	
		HardyCHROM™ MRSA	P14	Prepared Culture Media	N21	
		HardyCHROM™ VRE Agar	P19	Prepared Culture Media	N23	
		PDA with Lecithin and Tween®	P62	Prepared Culture Media	N37	
		R2A Agar with D/E	P54	Prepared Culture Media	N39	
		Sabdex Agar with Lecithin and Tween® 80	P25	Prepared Culture Media	N41	
		Template Sampling	SK701	Swabs & Transports	S12	
Tryptic Soy Agar with Lecithin and Tween® 80	P34	Prepared Culture Media	N48			
Fungi, Molds	Cultivation and Isolation	BHI Agar with Blood	X13	Prepared Culture Media	N4	
		BHI Agar with Chloramphenicol and Cycloheximide	X11	Prepared Culture Media	N5	
		BHI Agar Slant	L35	Prepared Culture Media	N4	
		BHI Agar with 5% Sheep Blood	W185	Prepared Culture Media	N4	
		BHI Agar with Blood and Gentamicin	X12	Prepared Culture Media	N4	
		Corn Meal Agar with Dextrose	W13	Prepared Culture Media	N13	
		Czapek Dox Agar Slant	L17	Prepared Culture Media	N14	
		Emmon's Sabdex Agar	W20	Prepared Culture Media	N40	
		Inhibitory Mold Agar with Chloramphenicol	X20	Prepared Culture Media	N23	
		Inhibitory Mold Agar with Gentamicin	W27	Prepared Culture Media	N24	
		Malt Extract Agar	W28	Prepared Culture Media	N28	
		Malt Extract Agar with 0.01% Chloramphenicol	W80	Prepared Culture Media	N28	
		m Green Yeast and Mold Agar	W63	Prepared Culture Media	N27	
		Mycobiotic Agar	L45, W50, X30	Prepared Culture Media	N32	
		Mycovue™ Slide Culture System	MV1	Prepared Culture Media	N32	
		Potato Dextrose Agar with TA	W96	Prepared Culture Media	N37	
		Potato Dextrose Agar	L90, W60	Prepared Culture Media	N37	
		Rose Bengal Agar with Chloramphenicol	H35, W87, U308	Prepared Culture Media	N40	
		Sabdex Agar	W70	Prepared Culture Media	N40	
		Sabdex Agar with Chloramphenicol and Tetracycline	W74	Prepared Culture Media	N41	
		SabHI Agar	X73	Prepared Culture Media	N41	
		Staining	BlueMount™	Z137	Stains & Reagents	Q9
	Lactophenol Cotton Blue		Z68	Rapid Tests & Reagents	P5	
	Potassium Hydroxide, 20% KOH		Z78	Stains & Reagents	Q9	
	Gardnerella vaginalis	Cultivation and Isolation	Gardnerella Selective Agar	G292	Prepared Culture Media	N19
			Starch Agar	G192	Prepared Culture Media	N44
			V Agar	A80	Prepared Culture Media	N51
V Agar with CNA			A81	Prepared Culture Media	N51	
V Agar/Starch Agar with BCP			J181	Prepared Culture Media	N51	
Identification		Hippurate Disks	DH4525	Rapid Tests & Reagents	P4	
		Hippurate Test with Ninhydrin	Z52	Rapid Tests & Reagents	P4	
		LAP Kit	DL10K	Rapid Tests & Reagents	P5	
		Metronidazole	DD8	Rapid Tests & Reagents	P7	
		Ninhydrin Reagent	Z70	Rapid Tests & Reagents	P8	
		Sulfonamide Disks	DD11	Rapid Tests & Reagents	P13	
Quality Control		<i>G. vaginalis</i> , ATCC® 14018™, KWIK-STIK™ DuoPak	0410P	Control Organisms	D9	
HACEK Group		Identification	CarboFerm™ II	Z99	Rapid Tests & Reagents	P3

Guide to Product Use



Organism	Usage	Product	Cat. no.	Section	Page
Giardia lamblia	Diagnostic Test	ImmunoCard STAT! [®] , Cryptosporidium/Giardia	750830	Direct Specimen Testing	H6
		Merifluor [®] Crypto/Giardia	250050	Parasitology	M15
		Para-Tect [™] Cryptosporidium/Giardia Direct Fluorescent Assay	MCCC675	Direct Specimen Testing	H6
		Para-Tect [™] Giardia lamblia Antigen Detection	MCCGA96	Direct Specimen Testing	H8
	Staining	Iodine, Dobell and O'Connor Trichrome Stain	Z66 T032	Rapid Tests & Reagents Stains & Reagents	P5 Q11
Quality Control	Control Slide, Intestinal Protozoa	SL1510	Parasitology	M14	
Haemophilus spp.	Cultivation and Isolation	Chocolate Agar	E14	Prepared Culture Media	N11
		Chocolate Agar with Bacitracin	E11	Prepared Culture Media	N12
		Chocolate Agar with Pyridoxal	E19	Prepared Culture Media	N12
		Chocolate Agar with Vancomycin	E16	Prepared Culture Media	N12
		Haemophilus ID Quadplate	J82	Prepared Culture Media	N20
		Peptic Digest Agar with Fildes & Nafcillin	G141	Prepared Culture Media	N35
		RTF Casman Selective Agar	A77	Prepared Culture Media	N40
	Identification	ALA Disks [™]	Z7081	Rapid Tests & Reagents	P1
		Beta Lactamase Test	Z51	Rapid Tests & Reagents	P1
		CarboFerm [™] II	Z99	Rapid Tests & Reagents	P3
		Haemophilus ID Quadplate for X and V Factors	J82	Prepared Culture Media	N20
		V Factor, HardyDisks [™]	Z7041	Rapid Tests & Reagents	P14
		X Factor, HardyDisk [™]	Z7031	Rapid Tests & Reagents	P15
	XV Factor, HardyDisk [™]	Z7051	Rapid Tests & Reagents	P15	
	Susceptibility	HTM Agar	G33, H07	Prepared Culture Media	N23
	Quality Control	<i>H. influenzae</i> , ATCC [®] 49247 [™] , KWIK-STIK [™] DuoPak	0647P	Control Organisms	D9
		<i>H. influenzae</i> , ATCC [®] 49766 [™] , KWIK-STIK [™] DuoPak	0919P	Control Organisms	D9
Helicobacter pylori	Identification	Indoxyl Acetate Disks	Z111	Rapid Tests & Reagents	P5
		ImmunoCard STAT! [®] HpSA [®]	750720	Direct Specimen Testing	H9
		QuickVue [®] One Step Helicobacter pylori gII [™]	W010	Direct Specimen Testing	H9
		Rapid Urea Medium	Z54	Prepared Culture Media	N51
		Urea Disks	DU7525	Rapid Tests & Reagents	P14
Lactobacillus spp.	Cultivation and Isolation	Lactobacilli MRS Agar	G117, G176	Prepared Culture Media	N24
		Lactobacilli MRS Agar with Cyclohexamide	G179, G214	Prepared Culture Media	N24
		Lactobacilli MRS Broth	K15, K17, K317	Prepared Culture Media	N24
	Quality Control	<i>L. acidophilus</i> , ATCC [®] 314 [™] , KWIK-STIK [™] DuoPak	0885P	Control Organisms	D11
Legionella spp.	Cultivation and Isolation	BCYE Agar	G07	Prepared Culture Media	N3
		BCYE Selective Agar with CAV Antibiotics	G08	Prepared Culture Media	N3
		BCYE Selective Agar with CCVC Antibiotics	G170	Prepared Culture Media	N3
		BCYE Selective Agar with GVPC Antibiotics	W169	Prepared Culture Media	N3
		BCYE with DGVP	G209	Prepared Culture Media	N3
	Environment	BCYE Agar	P251	Prepared Culture Media	N3
	Identification	Direct Fluorescent Monoclonal Antibody Test Kit	PL241	Direct Specimen Testing	H11
		Microgen [®] Legionella Latex Test	M45	Direct Specimen Testing	H11
		NOW [®] Legionella Urinary Antigen Test	852000	Direct Specimen Testing	H11
		Prolex [™] <i>L. pneumophila</i> Serogroup 1 Latex Kit	PL226	Direct Specimen Testing	H11
		SAS [™] Legionella Test	044020	Direct Specimen Testing	H11
	Quality Control	<i>L. pneumophila</i> , ATCC [®] 33152 [™] , KWIK-STIK [™] DuoPak	0211P	Control Organisms	D11

Guide to Product Use



Organism	Usage	Product	Cat. no.	Section	Page
Listeria spp.	Enrichment	Buffered Listeria Enrichment Broth Base	U233	Prepared Culture Media	N9
		Buffered Listeria Enrichment Broth Selective Supplement	K183	Prepared Culture Media	N9
		Fraser Broth	K98, U131	Prepared Culture Media	N19
		Fraser Broth Base	U57	Prepared Culture Media	N19
		Listeria Enrichment Broth, Modified	K57, U167	Prepared Culture Media	N9
	Cultivation and Isolation	UVM Broth	K96, U51	Prepared Culture Media	N51
		Oxford Medium, Modified	G46	Prepared Culture Media	N34
		HardyCHROM™ Listeria PALCAM	G317 G149	Prepared Culture Media Prepared Culture Media	N21 N34
	Identification	CAMP Spot Test	Z206	Rapid Tests & Reagents	P2
		Hippurate Test	Z52	Rapid Tests & Reagents	P4
		Microgen® Listeria ID Panel	MID67	Rapid Tests & Reagents	P7
		Microgen® Listeria Latex	F48	Rapid Tests & Reagents	P7
	Environment	ENVIROSWAB™	ENVSWB25	Swabs & Transports	S17
		InSite™ Rapid Environmental Listeria Test	IL100	Rapid Tests & Reagents	P5
		SDI-LIB Media for the detection of <i>Listeria</i> spp.	2500325	Swabs & Transports	S17
	Quality Control	<i>L. monocytogenes</i> , ATCC® 7644™, KWIK-STIK™ DuoPak	0398P	Control Organisms	D11
	Micrococcus spp.	Identification	Lysostaphin Differentiation Disks, HardyDisk™	Z112	Rapid Tests & Reagents
Quality Control		<i>M. luteus</i> , ATCC® 4698™, KWIK-STIK™ DuoPak	0242P	Control Organisms	D12
Moraxella catarrhalis	Cultivation and Isolation	Moraxella catarrhalis Selective Agar	A135	Prepared Culture Media	N30
	Identification	CarboFerm™ I	Z98	Rapid Tests & Reagents	P2
		CatScreen™ Rapid Test	Z110	Rapid Tests & Reagents	P3
Quality Control	<i>M. catarrhalis</i> , ATCC® 25238™, KWIK-STIK™ DuoPak	0951P	Control Organisms	D12	
	<i>M. catarrhalis</i> , ATCC® 25240™, KWIK-STIK™ DuoPak	0403P	Control Organisms	D12	
Mycobacterium spp.	Digestion and Decontamination	MucoGest™ 50	Z130	Rapid Tests & Reagents	P8
		NALC Powder	Z60	Rapid Tests & Reagents	P8
		TB Prep Kit Red™	Z158	Rapid Tests & Reagents	P14
		TB Prep Kit™	Z160	Rapid Tests & Reagents	P14
	Cultivation and Isolation	LJ Deep	C27	Prepared Culture Media	N26
		LJ Gruft	C23, X23	Prepared Culture Media	N26
		LJ Media	X22	Prepared Culture Media	N26
		LJ Selective Slant	C25, X34	Prepared Culture Media	N26
		LJ Slant	C21	Prepared Culture Media	N26
		LJ without Glycerol (with Pyruvate) Slant	C31	Prepared Culture Media	N26
		LJ with Pyruvate	C28, X19	Prepared Culture Media	N26
		Middlebrook 7H10 Agar	C34, W30, X26	Prepared Culture Media	N29
		Middlebrook 7H11 Agar	J75, C36, W35	Prepared Culture Media	N29
	Middlebrook 7H11 Selective	X28	Prepared Culture Media	N30	
	Identification	AFB Nitrate Reagent A	Z124	Rapid Tests & Reagents	P1
		AFB Nitrate Reagent B	Z125	Rapid Tests & Reagents	P1
		Arylsulfatase Broth	K93, K94	Prepared Culture Media	N2
		LJ with Sodium Chloride Slant	C29	Prepared Culture Media	N26
		Middlebrook 7H9 Broth with Tween® 80	C62	Prepared Culture Media	N30
		Pyrazinamidase Agar	K158	Prepared Culture Media	N39
		Urea Disks	DU7525	Rapid Tests & Reagents	P14
	Staining	Acid Fast Fluorescent Stain Set M Auramine	418K	Stains & Reagents	Q3
		Acid Fast Stain Kit, Kinyoun	AF900	Stains & Reagents	Q1
		Acid Fast Stain Kit, Ziehl-Neelsen	484K	Stains & Reagents	Q2
		Bovine Serum Albumin 0.2%	Z81	Rapid Tests & Reagents	P2
		Methylene Blue	MB008	Featured Products	T27
	Quality Control	Q-Slide™ AFB (Acid Fast Bacillus)	Z301	Control Organisms	D42
		<i>M. tuberculosis</i> , ATCC® 25177™, KWIK-STIK™ DuoPak	0112P	Control Organisms	D13

Guide to Product Use



Organism	Usage	Product	Cat. no.	Section	Page
Mycoplasma spp.	Cultivation and Isolation	A8 Agar	G02	Prepared Culture Media	N1
		Mycoplasma Agar (bovine mastitis)	G102	Prepared Culture Media	N32
		Mycoplasma Agar with Cefoperazone (bovine mastitis)	G212	Prepared Culture Media	N32
		Mycoplasma Broth (bovine mastitis)	R102	Prepared Culture Media	N32
		PPLO Selective Agar	G04, G105	Prepared Culture Media	N38
		SP4 Agar with Arginine	G32	Prepared Culture Media	N43
		SP4 Agar with Glucose	G21	Prepared Culture Media	N43
		SP4 Broth Plain	R84	Prepared Culture Media	N43
	Identification	ImmunoCard®, Mycoplasma	709030	Direct Specimen Testing	H13
	Transport	CVM Flocked Swab Kit Female	R96KF	Direct Specimen Testing	H2
		CVM Flocked Swab Kit Male	R96KM	Direct Specimen Testing	H2
		CVM Flocked Swab Kit Nasopharyngeal	R96KNP	Direct Specimen Testing	H2
		CVM Transport Media	R96	Direct Specimen Testing	H2
Quality Control	<i>U. parvum</i> , ATCC® 27813™, KWIK-STIK™ DuoPak	0151P	Control Organisms	D18	
Naegleria spp.	Cultivation and Isolation	Non-Nutrient Agar	G225	Prepared Culture Media	N33
		Page's Saline	R225	Prepared Culture Media	N34
	Maintenance of Stock Cultures	Naegleria Broth	K325	Prepared Culture Media	N33
	Staining	Modified Field's Stain Kit	Z15	Stains & Reagents	Q16
Neisseria gonorrhoeae	Cultivation and Isolation	Chocolate/Martin Lewis with Lincomycin	J44	Prepared Culture Media	N11
		Chocolate/Modified Thayer Martin	J72	Prepared Culture Media	N11
		Jembec Martin Lewis	1737	Prepared Culture Media	N29
		Martin Lewis with Lincomycin	E39	Prepared Culture Media	N29
		Martin Lewis with Lincomycin with Pill Pocket	E31	Prepared Culture Media	N29
		CO ₂ Tablet (for Martin Lewis with Lincomycin)	CO2TAB	Anaerobic Microbiology	A14
		Thayer Martin Modified Agar	E30	Prepared Culture Media	N46
		Transgrow	X52	Prepared Culture Media	N47
	Identification	Beta Lactamase Test	Z51	Rapid Tests & Reagents	P1
		CarboFerm™ I Neisseria	Z98	Rapid Tests & Reagents	P2
		GonoGen II	89102024, 89102040	See Web Catalog	—
		OxiDrops™	Z119	Rapid Tests & Reagents	P9
		OxiStrips™	Z93	Rapid Tests & Reagents	P10
		WeeTab Neisseria	K9187	Rapid Tests & Reagents	P15
	Quality Control	<i>N. gonorrhoeae</i> , ATCC® 43069™, KWIK-STIK™ DuoPak	0429P	Control Organisms	D13
		<i>N. gonorrhoeae</i> , ATCC® 49226™, KWIK-STIK™ DuoPak	0648P	Control Organisms	D13
	Nocardia spp.	Cultivation and Isolation	BCYE Agar	G07	Prepared Culture Media
Blood Agar			A10	Prepared Culture Media	N6
Chocolate Agar			E14	Prepared Culture Media	N11
Thayer Martin Modified Agar with VCNT			E30	Prepared Culture Media	N46
Identification		Nocardia Hydrolysis Kit	Z604	Prepared Culture Media	N33
Quality Control		<i>N. brasiliensis</i> , ATCC® 19296™, KWIK-STIK™ DuoPak	0866P	Control Organisms	D13
Parasitology	Quality Control	Control Slide, Cryptosporidium	SL4510	Control Organisms	D43
		Control Slide, Microsporidium	SL8010	Control Organisms	D43
		Control Slide, SAF-Preserved Fecal Smear Control	SL1510	Control Organisms	D43
		Control Slide, Zinc PVA-Preserved Fecal Smear Control	SL1010	Control Organisms	D43
	Identification	ParaCon™	PK48	Parasitology	M5
	Staining	Trichrome Stain	T032	Stains & Reagents	Q2
Trichrome Starter Kit		VXF019	Parasitology	M8	
Plasmodium spp. (malaria)	Quality Control	Control Slide, Blood Parasite	SL9010	Control Organisms	D43
	Staining	Giemsa Blood Stain	591A16	Stains & Reagents	Q16
Stat Stain, Rapid Wright-Giemsa		SS016	Stains & Reagents	Q16	
Pneumocystis	Diagnostic Test	Merifluor®, Pneumocystis	222030	Direct Specimen Testing	H14
	Quality Control	Control Slide, Blood Parasites	SL9010	Control Organisms	D43

Guide to Product Use



Organism	Usage	Product	Cat. no.	Section	Page
<i>Pseudomonas aeruginosa</i>	Cultivation and Isolation	Cetrimide Select Agar (CET)	G18, L34	Prepared Culture Media	N11
		CET/MAC/VJ Agar	J314	Prepared Culture Media	N11
		HardyCHROM™ Urine Biplate	J100	Prepared Culture Media	N21
		HardyCHROM™ UTI	G313	Prepared Culture Media	N22
		Pseudomonas Isolation Agar	G219	Prepared Culture Media	N38
	Identification	OxiStrips™	Z93	Rapid Tests & Reagents	P10
	Water	Asparagine Broth	K170, K173	Prepared Culture Media	N2
		m PA Agar	G133	Prepared Culture Media	N27
		m PA-C Agar	G150	Prepared Culture Media	N27
	Quality Control	<i>P. aeruginosa</i> , ATCC® 27853™, KWIK-STIK™ DuoPak	0353P	Control Organisms	D44
<i>Salmonella</i> spp.	Enrichment	Buffered Peptone Water	K195, U142	Prepared Culture Media	N9
		GN Broth	K01	Prepared Culture Media	N19
		Hajna Tetrathionate (TT) Broth	K289	Prepared Culture Media	N20
		Iodine-Iodide Solution	Z129, Z139	Prepared Culture Media	N24
		Selenite Cystine Broth	K69	Prepared Culture Media	N42
		Tetrathionate Broth	K65, K165, U165	Prepared Culture Media	N46
		Tetrathionate Broth with Brilliant Green	K164	Prepared Culture Media	N46
		Tetrathionate Broth with Brilliant Green Bile	K273	Prepared Culture Media	N46
		Cultivation and Isolation	BG Sulfa Agar	G87	Prepared Culture Media
	Brilliant Green Agar		G75	Prepared Culture Media	N7
	HardyCHROM™ Salmonella		G309	Prepared Culture Media	N21
	HardyCHROM™ SS Agar		G327	Prepared Culture Media	N22
	HE Agar		G63	Prepared Culture Media	N23
	HE Agar with Novobiocin		G242	Prepared Culture Media	N23
	HE Agar/SS Agar		J139	Prepared Culture Media	N23
	LIA Agar, Modified		G86	Prepared Culture Media	N25
	MacConkey Agar		G35	Prepared Culture Media	N28
	SS Agar		G50	Prepared Culture Media	N43
	Tetrathionate Broth		K165	Prepared Culture Media	N46
	XLD Agar		G65	Prepared Culture Media	N52
	XLT4 Agar		G165, G277	Prepared Culture Media	N53
	Environment		EnviroBootie™, EnviroBootie™ II	EB20, EB20M, EB60K	Swabs & Transports
	Identification	Antisera	MULTIPLE	Antisera	————
		BioLine Salmonella ELISA Test Kit	CHB0081	See Web Catalog	————
		EnteroScreen 4™	L225	Prepared Culture Media	N16
		LIA Slant	L25	Prepared Culture Media	N26
		Microgen® GN- ID Panel A	MID64	Rapid Tests & Reagents	P7
		Microgen® Salmonella Latex Test	M42	Rapid Tests & Reagents	P12
		ONPG-PDA Test	R93	Rapid Tests & Reagents	P9
		Triple Sugar Iron (TSI) Agar Slant	L16, R32	Prepared Culture Media	N47
	Quality Control	<i>S. choleraesuis</i> , ATCC® 10708™, KWIK-STIK™ DuoPak	0902P	Control Organisms	D15
		<i>S. enteritidis</i> , ATCC® 13076™ KWIK-STIK™ DuoPak	0345P	Control Organisms	D15
		<i>S. enterica</i> , ATCC® 14028™, EZ-Accu Shot	0363A	Control Organisms	D15
<i>Shigella</i> spp.	Enrichment	GN Broth	K01	Prepared Culture Media	N19
	Cultivation and Isolation	HardyCHROM™ SS Agar	G327	Prepared Culture Media	N22
		HE Agar	G63	Prepared Culture Media	N23
		MacConkey Agar	G35	Prepared Culture Media	N28
		SS Agar	G50	Prepared Culture Media	N43
	Identification	Antisera	MULTIPLE	Antisera	————
		EnteroScreen 4™	L225	Prepared Culture Media	N16
		LIA Slant	L25	Prepared Culture Media	N26
		Microgen® GN- ID Panel A	MID64	Rapid Tests & Reagents	P7
	Triple Sugar Iron (TSI) Agar Slant	L16, R32	Prepared Culture Media	N47	
	Quality Control	<i>S. flexneri</i> , ATCC® 12022™, KWIK-STIK™ DuoPak	0356P	Control Organisms	D16

Guide to Product Use



Organism	Usage	Product	Cat. no.	Section	Page
<i>Staphylococcus aureus</i>	Cultivation and Isolation	Baird-Parker Agar	G96	Prepared Culture Media	N2
		BEA with Azide/CNA Agar	J66	Prepared Culture Media	N4
		Blood Agar, 5%	A10	Prepared Culture Media	N6
		Columbia CNA	A50	Prepared Culture Media	N12
		HardyCHROM™ <i>Staphylococcus aureus</i>	G311	Prepared Culture Media	N22
		Mannitol Salt Agar	G40	Prepared Culture Media	N31
		Rose Agar	A66	Prepared Culture Media	N40
		Staph 110 Agar	G52	Prepared Culture Media	N44
	Vogel and Johnson Agar	G193	Prepared Culture Media	N52	
	Identification	Catalase Reagent	Z62	Rapid Tests & Reagents	P2
		Coagulase Cryo™	Z202	Rapid Tests & Reagents	P3
		CoaguStaph™	Z003, Z020	Rapid Tests & Reagents	P3
		DNase Agar with Toluidine Blue	G24	Prepared Culture Media	N15
		StaphTEX™ Blue, Latex Agglutination Test	ST50, ST200	Rapid Tests & Reagents	P13
	Water	m <i>Staphylococcus</i> Broth	K213	Prepared Culture Media	N27
Quality Control	Control Slide, <i>S. aureus</i>	SL0110	Control Organisms	D43	
	<i>S. aureus</i> , ATCC® 6538™, KWIK-STIK™ DuoPak	0485P	Control Organisms	D44	
	<i>S. aureus</i> , ATCC® 29213™, KWIK-STIK™ DuoPak	0365P	Control Organisms	D16	
<i>Staphylococcus aureus</i>, Methicillin Resistant	Cultivation and Isolation	HardyCHROM™ MRSA	G307	Prepared Culture Media	N21
		MSA with Oxacillin	G97	Prepared Culture Media	N31
		MSA/MSA with Oxacillin	J197	Prepared Culture Media	N31
	Identification	Alere™ PBP2a	892000	Rapid Tests & Reagents	P1
		MRSA Latex	DR900A	Rapid Tests & Reagents	P8
	Susceptibility	MRSA Screen Plate	G47	Prepared Culture Media	N31
	Environment	HardyCHROM™ MRSA	P14	Prepared Culture Media	N21
	Quality Control	<i>S. aureus</i> , ATCC® 25923™, KWIK-STIK™ DuoPak	0360MSSA	Control Organisms	D19
		<i>S. aureus</i> , ATCC® 43300™, KWIK-STIK™ DuoPak	0852P	Control Organisms	D16
	<i>Staphylococcus lugdunensis</i>	Identification	PYR Test Kit	Z75	Rapid Tests & Reagents
Rapid Ornithine			K279	Rapid Tests & Reagents	P11
<i>Streptococcus</i> group A	Cultivation and Isolation	Group A Beta Strep Agar	A72	Prepared Culture Media	N20
		Selective Beta Strep Agar/Selective Beta Strep Agar	J73	Prepared Culture Media	N20
	Identification	Bacitracin 0.04 units, HardyDisks™	Z7021	Rapid Tests & Reagents	P1
		PYR Test Kit	Z75	Rapid Tests & Reagents	P10
		QuickVue+ Strep A Test	Z0122	See Web Catalog	
		StrepQuick™	Z122	Rapid Tests & Reagents	P13
	Quality Control	StrepPRO™ Streptococcal Grouping Kit	PL030HD	Rapid Tests & Reagents	P13
		<i>S. agalactiae</i> , ATCC® 12386™, KWIK-STIK™ DuoPak	0439P	Control Organisms	D17
	<i>S. pyogenes</i> , ATCC® 19615™, KWIK-STIK™ DuoPak	0385P	Control Organisms	D18	
	<i>Streptococcus</i> group B	Cultivation and Isolation	CAMP Spot Test Reagent	Z206	Rapid Tests & Reagents
GBS Detect™			A300	Prepared Culture Media	N19
Granada Medium			G123	Prepared Culture Media	N20
LIM Broth			L57	Prepared Culture Media	N26
Selective Strep Agar			A70	Prepared Culture Media	N42
Strep B Carrot Broth™			Z140	Prepared Culture Media	N44
Identification		StrepPRO™ Streptococcal Grouping Kit	PL030HD	Rapid Tests & Reagents	P13
Quality Control		<i>S. agalactiae</i> , ATCC® 13813™, KWIK-STIK™ DuoPak	0370P	Control Organisms	D17
		<i>S. agalactiae</i> , ATCC® 12386™, KWIK-STIK™ DuoPak	0439P	Control Organisms	D17
<i>Streptococcus anginosus</i>		Identification	Rapid Anginosus ID Kit	Z14	Rapid Tests & Reagents
	Quality Control	<i>S. anginosus</i> , NCTC 10713, KWIK-STIK™ DuoPak	0389P	Control Organisms	D17
<i>Streptococcus</i> (Beta-hemolytic) A, B, C, D, F, G	Identification	StrepPRO™ Streptococcal Grouping Kit	PL030HD	Rapid Tests & Reagents	P13

Guide to Product Use



Organism	Usage	Product	Cat. no.	Section	Page	
<i>Streptococcus pneumoniae</i>	Cultivation and Isolation	Blood Agar Plate, 5%	A10	Prepared Culture Media	N6	
		Selective Strep Agar	A70	Prepared Culture Media	N42	
		Streptococcus pneumoniae Selective Agar	A133	Prepared Culture Media	N45	
	Identification	Optochin, HardyDisk™	Z7011	Rapid Tests & Reagents	P9	
		Bile Spot Reagent	Z61	Rapid Tests & Reagents	P2	
Quality Control	<i>S. pneumoniae</i> , ATCC® 6305™, KWIK-STIK™ DuoPak	0447P	Control Organisms	D18		
	<i>S. pneumoniae</i> , ATCC® 49619™, KWIK-STIK™ DuoPak	0947P	Control Organisms	D18		
<i>Toxoplasmosis spp.</i>	Diagnostic Test	Mast Toxoreagent™	RST701	Direct Specimen Testing	H19	
		Toxoplasma Test, TPM®	40D3	Direct Specimen Testing	H19	
<i>Treponema pallidum (syphilis)</i>	Diagnostic Test	RPR Kit	300625550	Direct Specimen Testing	H18	
		Serodia-Tp-Pa Agglutination Test	1626	Direct Specimen Testing	H18	
		VDRL Antigen Kit	950010	Direct Specimen Testing	H18	
<i>Trichomonas spp.</i>	Cultivation and Isolation	Diamonds Medium, Modified	K02	Prepared Culture Media	N15	
		InPouch™ TV	MULTIPLE	Parasitology	M14	
<i>Ureaplasma</i>	Cultivation and Isolation	A8 Agar	G02	Prepared Culture Media	N1	
<i>Vibrio spp.</i>	Cultivation and Isolation	Alkaline Peptone Water	K03	Prepared Culture Media	N1	
		TCBS Agar	Q52, G55	Prepared Culture Media	N45	
	Identification	O129 Disk	DD14, DD15	Rapid Tests & Reagents	P9	
	Vibrio Cholerae Latex Kit	296092	See Web Catalog	—		
Quality Control	<i>V. parahaemolyticus</i> , ATCC® 17802™, KWIK-STIK™ DuoPak	0818P	Control Organisms	D18		
<i>Yersinia enterocolitica</i>	Cultivation and Isolation	CIN Agar	G20	Prepared Culture Media	N12	
		CIN Agar/MacConkey Agar with Sorbitol	J49	Prepared Culture Media	N12	
		Microgen® GN-ID	MID641	Rapid Tests & Reagents	P7	
		Yersinia Selective Agar Base, CRITERION™	C5391	Dehydrated Culture Media	G9	
Virus	Diagnostic Test	Adeno virus Type 1	301215	Antisera	B14	
		Biokit HSV-2 Rapid Test Kit	300028001	Direct Specimen Testing	H9	
		Color Mono Kit	30001005	Direct Specimen Testing	H12	
		ImmunoCard STAT!® Rotavirus	750030	Direct Specimen Testing	H16	
		Merifluor, VZV	500113	Direct Specimen Testing	H19	
		QuickVue® Influenza A+B	20183	Direct Specimen Testing	H9	
		NOW® RSV	430000	Direct Specimen Testing	H16	
		Rotascreen® Latex Test	M80	Direct Specimen Testing	H16	
		Rubajet Rubella Latex Agglutination Slide Kit	300024011	Direct Specimen Testing	H17	
		SAS™ Rota Test	069020	Direct Specimen Testing	H16	
		SAS™ RSV Alert	041130	Direct Specimen Testing	H16	
		Transport	CVM Flocked Swab Kit Female	R96KF	Direct Specimen Testing	H2
			CVM Flocked Swab Kit Male	R96KM	Direct Specimen Testing	H2
	CVM Flocked Swab Kit Nasopharyngeal		R96KNP	Direct Specimen Testing	H9	
	CVM Transport Media		R96	Direct Specimen Testing	H9	
	Miscellaneous	Stock Culture Preservation	CryoSavers™ Cryogenic Vials	MULTIPLE	Prepared Cultured Media	N54





800.266.2222
805.346.2760 (fax)
www.HardyDiagnostics.com
Sales@HardyDiagnostics.com



The Quality Management System at Hardy Diagnostics manufacturing facility is certified to ISO 13485.

ISO 13485
FM 572526

Santa Maria, CA | Centralia, WA | Salt Lake City, UT | Phoenix, AZ | Dallas, TX | Springboro, OH | Lake City, FL | Albany, NY