



SUSCEPTIBILITY



TESTING



# SUSCEPTIBILITY TESTING

## Table of Contents

- R1 HardyDisks™ AST
- R6 HardyDisks™ ID
- R7 BBL™ Sensi-Disc™
- R10 Accessory Products for AST
- R11 Culture Media for AST



### Headquarters

1430 West McCoy Lane  
 Santa Maria, CA 93455  
 800.266.2222 : phone  
 805.346.2760 : fax  
 Sales@HardyDiagnostics.com  
 www.HardyDiagnostics.com

### Distribution Centers

Santa Maria, California  
 Centralia, Washington  
 Salt Lake City, Utah  
 Phoenix, Arizona  
 Dallas, Texas  
 Springboro, Ohio  
 Lake City, Florida  
 Albany, New York



ISO 13485  
 FM 572526

The Quality Management System at the Hardy Diagnostics manufacturing facility is certified to ISO 13485.

Copyright © 2011 Hardy Diagnostics

# HardyDisks™ AST

*Now you have a competitive choice!*

Hardy Diagnostics offers an extensive selection of antibiotic sensitivity disks for the Kirby-Bauer disk diffusion test.



Mueller Hinton Agar, Cat. no. G45



## Features & Benefits

- ⊗ A complete line of all commonly used antibiotics.
- ⊗ Compatible with the BBL™ Sensi-Disc™ disk dispensers (6, 8, or 12-place).
- ⊗ Spring-loaded in cartridges of 50 disks for smooth and reliable dispensing.
- ⊗ Can also be used with the metal single disk dispenser (Cat. no. 260457).
- ⊗ Available in single cartridges or packs of five.
- ⊗ Improved to prevent jamming and tearing.




## Innovative Packaging!

HardyDisks™ AST cartridges are packaged in a new canister that features a desiccant built into the lid.

## New HardyDisk Last Disk Recognition!

The last disk in each cartridge is clearly marked with dark lines to signify that a new cartridge needs to be loaded into the dispenser. Prevents you from dispensing "missing disks"



## Compatible with Common Dispensers

The HardyDisks™ AST cartridges are compatible with the BBL™ single or multiple disk dispensers. Each cartridge contains 50 disks that are spring-loaded for smooth and reliable dispensing.

# HardyDisks™ AST



Product Name	Code	Unit of Measure	Cat. no.
Amikacin	AN-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z8011
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z8015
Amoxicillin/Clavulanic Acid, Augmentin	AmC-30	20µg/10µg, 50 disks/cartridge, 1 cartridge/canister	Z8021
		20µg/10µg, 50 disks/cartridge, 5 cartridges/canister	Z8025
Ampicillin/Sulbactam	SAM-20	10µg/10µg, 50 disks/cartridge, 1 cartridge/canister	Z8041
		10µg/10µg, 50 disks/cartridge, 5 cartridges/canister	Z8045
Ampicillin	AM-10	10µg, 50 disks/cartridge, 1 cartridge/canister	Z8031
		10µg, 50 disks/cartridge, 5 cartridges/canister	Z8035
Azithromycin	AZM-15	15µg, 50 disks/cartridge, 1 cartridge/canister	Z8051
		15µg, 50 disks/cartridge, 5 cartridges/canister	Z8055
Bacitracin	B-10	10µg, 50 disks/cartridge, 1 cartridge/canister	Z9271
		10µg, 50 disks/cartridge, 5 cartridges/canister	Z9275
Carbenicillin	CB-100	100µg, 50 disks/cartridge, 1 cartridge/canister	Z8111
		100µg, 50 disks/cartridge, 5 cartridges/canister	Z8115
Cefaclor	CEC-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z8121
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z8125
Cefazolin	CZ-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z8141
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z8145
Cefixime	CFM-5	5µg, 50 disks/cartridge, 1 cartridge/canister	Z8181
		5µg, 50 disks/cartridge, 5 cartridges/canister	Z8185
Cefotaxime	CTX-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z8221
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z8225
Cefotetan	CTT-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z8231
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z8235

Product Name	Code	Unit of Measure	Cat. no.
Cefoxitin	FOX-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z8241
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z8245
Ceftaroline	CPT-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z9391
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z9395
Ceftazidime	CAZ-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z8271
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z8275
Ceftriaxone	CRO-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z8311
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z8315
Cefuroxime	CXM-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z8321
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z8325
Cephalothin	CF-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z8331
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z8335
Chloramphenicol	C-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z8341
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z8345
Cinoxacin	CIN-100	100µg, 50 disks/cartridge, 1 cartridge/canister	Z8351
		100µg, 50 disks/cartridge, 5 cartridges/canister	Z8355
Ciprofloxacin	CIP-5	5µg, 50 disks/cartridge, 1 cartridge/canister	Z8361
		5µg, 50 disks/cartridge, 5 cartridges/canister	Z8365
Clarithromycin	CLR-15	15µg, 50 disks/cartridge, 1 cartridge/canister	Z8371
		15µg, 50 disks/cartridge, 5 cartridges/canister	Z8375
Clindamycin	CC-2	2µg, 50 disks/cartridge, 1 cartridge/canister	Z8391
		2µg, 50 disks/cartridge, 5 cartridges/canister	Z8395
Colistin Sulfate	CT	10µg, 50 disks/cartridge, 1 cartridge/canister	Z8411
Doripenem	DOR-10	30µg, 50 disks/cartridge, 1 cartridge/canister	Z9371
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z9375
Doxycycline	D-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z8441
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z8445
Erythromycin	E-15	15µg, 50 disks/cartridge, 1 cartridge/canister	Z8471
		15µg, 50 disks/cartridge, 5 cartridges/canister	Z8475





Product Name	Code	Unit of Measure	Cat. no.
Gentamicin	GM-10	10µg, 50 disks/cartridge, 1 cartridge/canister	Z8541
		10µg, 50 disks/cartridge, 5 cartridges/canister	Z8545
Imipenem	IPM-10	10µg, 50 disks/cartridge, 1 cartridge/canister	Z8571
		10µg, 50 disks/cartridge, 5 cartridges/canister	Z8575
Kanamycin	K-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z8601
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z8605
Levofloxacin	LVX-5	5µg, 50 disks/cartridge, 1 cartridge/canister	Z8611
		5µg, 50 disks/cartridge, 5 cartridges/canister	Z8615
Metallo-β-Lactamase: Set of 2 Cartridges			
Imipenem	IMI10	10µg, 50 disks/cartridge, 1 cartridge/canister	Z155
Imipenem + EDTA	IM10E	10µg/750µg, 50 disks/cartridge, 1 cartridge/canister Not compatible with the BBL Single Disk Dispensers	
Nalidixic Acid	NA-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z8731
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z8735
Nitrofurantoin	Fm-300	300µg, 50 disks/cartridge, 1 cartridge/canister	Z8791
		300µg, 50 disks/cartridge, 5 cartridges/canister	Z8795
Norfloxacin	NOR-10	10µg, 50 disks/cartridge, 1 cartridge/canister	Z8811
		10µg, 50 disks/cartridge, 5 cartridges/canister	Z8815
Novobiocin	NB-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z8831
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z8835
Ofloxacin	OFX-5	5µg, 50 disks/cartridge, 1 cartridge/canister	Z8841
		5µg, 50 disks/cartridge, 5 cartridges/canister	Z8845
Oxacillin	OX-1	1µg, 50 disks/cartridge, 1 cartridge/canister	Z8851
		1µg, 50 disks/cartridge, 5 cartridges/canister	Z8855
Penicillin	P-10	10µg, 50 disks/cartridge, 1 cartridge/canister	Z8881
		10µg, 50 disks/cartridge, 5 cartridges/canister	Z8885
Piperacillin	PIP-100	100µg, 50 disks/cartridge, 1 cartridge/canister	Z8901
		100µg, 50 disks/cartridge, 5 cartridges/canister	Z8905

Product Name	Code	Unit of Measure	Cat. no.
Piperacillin/Tazobactam	TZP-110	100µg/10µg, 50 disks/cartridge, 1 cartridge/canister	Z8911
		100µg/10µg, 50 disks/cartridge, 5 cartridges/canister	Z8915
Streptomycin	S-10	10µg, 50 disks/cartridge, 1 cartridge/canister	Z8991
		10µg, 50 disks/cartridge, 5 cartridges/canister	Z8995
Sulfamethoxazole/Trimethoprim	SXT	1.25µg/23.75µg, 50 disks/cartridge, 1 cartridge/canister	Z9011
		1.25µg/23.75µg, 50 disks/cartridge, 5 cartridges/canister	Z9015
Sulfisoxazole	G-0.25	250µg, 50 disks/cartridge, 1 cartridge/canister	Z9071
		250µg, 50 disks/cartridge, 5 cartridges/canister	Z9075
Telavancin	TLV-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z9381
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z9385
Tetracycline	Te-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z9121
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z9125
Ticarcillin	TIC-75	75µg, 50 disks/cartridge, 1 cartridge/canister	Z9141
		75µg, 50 disks/cartridge, 5 cartridges/canister	Z9145
Tigecycline	TGC-15	15µg, 50 disks/cartridge, 1 cartridge/canister	Z9361
Timentin, (Ticarcillin with Clavulanic Acid)	TIM-85	75µg/10µg, 50 disks/cartridge, 1 cartridge/canister	Z9161
		75µg/10µg, 50 disks/cartridge, 5 cartridges/canister	Z9165
Tobramycin	NN-10	10µg, 50 disks/cartridge, 1 cartridge/canister	Z9171
		10µg, 50 disks/cartridge, 5 cartridges/canister	Z9175
Trimethoprim	TMP-5	5µg, 50 disks/cartridge, 1 cartridge/canister	Z9181
		5µg, 50 disks/cartridge, 5 cartridges/canister	Z9185
Vancomycin	Va-30	30µg, 50 disks/cartridge, 1 cartridge/canister	Z9241
		30µg, 50 disks/cartridge, 5 cartridges/canister	Z9245

# HardyDisks™ ID

Product Name	Unit of Measure	Cat. no.
ALA Disks™	50 disks/cartridge, 1 cartridge	Z7081
Bacitracin, 0.04 unit	50 disks/cartridge, 1 cartridge	Z7021
	50 disks/cartridge, 5 cartridges	Z7025
Bile Disk	50 disks/cartridge, 1 cartridge	Z7091
Blank Disk, 0.25in	50 disks/cartridge, 1 cartridge	Z7121
Kanamycin, 1000µg	50 disks/cartridge, 1 cartridge	Z7191
Lysostaphin	50 disks/vial	Z112
Nitrate Disks	50 disks/cartridge, 1 cartridge	Z7071
Nitrocef Disks	50 disks/vial	Z7301
Novobiocin, 5µg, NB-5	50 disks/cartridge, 1 cartridge	Z7291
ONPG Disks	50 disks/cartridge, 1 cartridge	Z7321
Optochin	50 disks/cartridge, 1 cartridge	Z7011
	50 disks/cartridge, 5 cartridges	Z7015
SPS Disks, 1mg	50 disks/cartridge, 1 cartridge	Z7381
Urea/PDA Disks	50 disks/cartridge, 1 cartridge	Z7331
Vancomycin, 5µg, Va-5	50 disks/cartridge, 1 cartridge	Z7501
V Factor Disks	50 disks/cartridge, 1 cartridge	Z7041
X Factor Disks	50 disks/cartridge, 1 cartridge	Z7031
XV Factor Disks	50 disks/cartridge, 1 cartridge	Z7051



# BBL™ Sensi-Disc™ by BD

Product Name	Code	Unit of Measure	Cat. no.
Amdinocillin	AMD-10	10µg, 50 disks, 1 cartridge	231638
Amikacin	AN-30	30µg, 50 disks, 1 cartridge	231596
Ampicillin	AM-2	2µg, 50 disks, 1 cartridge	230701
Ampicillin	AM-10	10µg, 50 disks, 1 cartridge	230705
Augmentin	AmC-30	30µg, 50 disks, 1 cartridge	231628
Azithromycin	AZM-15	15µg, 50 disks, 1 cartridge	231681
Azlocillin	AZ-75	75µg, 50 disks, 1 cartridge	231624
Aztreonam	ATM-30	30µg, 50 disks, 1 cartridge	231640
Bacitracin	B-2	2µg, 50 disks, 1 cartridge	230717
Bacitracin	B-10	10µg, 50 disks, 1 cartridge	230721
Carbenicillin	CB-100	100µg, 50 disks, 1 cartridge	231235
Cefaclor	CEC-30	30µg, 50 disks, 1 cartridge	231652
Cefamandole	MA-30	30µg, 50 disks, 1 cartridge	231594
Cefazolin	CZ-30	30µg, 50 disks, 1 cartridge	231592
Cefepime	FEP-30	30µg, 50 disks, 1 cartridge	231695
Cefixime	CFM-5	5µg, 50 disks, 1 cartridge	231663
Cefmetazole	CMZ-30	30µg, 50 disks, 1 cartridge	231665
Cefonicid	CID-30	30µg, 50 disks, 1 cartridge	231642
Cefoperazone	CFP-75	75µg, 50 disks, 1 cartridge	231612
Cefotaxime	CTX-30	30µg, 50 disks, 1 cartridge	231606
Cefotaxime/Clavulanic acid	CTX/CLA	30µg/10µg, 50 disks, 1 cartridge	231751
Cefotetan	CTT-30	30µg, 50 disks, 1 cartridge	231655
Cefoxitin	FOX-30	30µg, 50 disks, 1 cartridge	231590
Cefpodoxime	CPD-10	10µg, 50 disks, 1 cartridge	231673
Cefprozil	CPR-30	30µg, 50 disks, 1 cartridge	231683
Ceftazidime	CAZ-30	30µg, 50 disks, 1 cartridge	231632
Ceftazidime/Clavulanic acid	CAZ/CLA	30µg/10µg, 50 disks, 1 cartridge	231753
Ceftizoxime	ZOX-30	30µg, 50 disks, 1 cartridge	231622
Ceftriaxone	CR0-30	30µg, 50 disks, 1 cartridge	231634
Cefuroxime	CXM-30	30µg, 50 disks, 1 cartridge	231620
Cephalothin	CF-30	30µg, 50 disks, 1 cartridge	230725
Chloramphenicol	C-5	5µg, 50 disks, 1 cartridge	230729
Chloramphenicol	C-30	30µg, 50 disks, 1 cartridge	230733
Cinoxacin	CIN-100	100µg, 50 disks, 1 cartridge	231598
Ciprofloxacin	CIP-5	5µg, 50 disks, 1 cartridge	231657
Clarithromycin	CLR-15	15µg, 50 disks, 1 cartridge	231677

Product Name	Code	Unit of Measure	Cat. no.
Clindamycin	CC-2	2µg, 50 disks, 1 cartridge	231213
Cloxacillin	CX-1	1µg, 50 disks, 1 cartridge	230737
Colistin	CL-10	10µg, 50 disks, 1 cartridge	230749
Doxycycline	D-30	30µg, 50 disks, 1 cartridge	230777
Enoxacin	ENX-10	10µg, 50 disks, 1 cartridge	231661
Enrofloxacin	ENO-5	5µg, 50 disks, 1 cartridge	231717
Ertapenem	ETP-10	10µg, 50 disks, 1 cartridge	232174
Erythromycin	E-2	2µg, 50 disks, 1 cartridge	230789
Erythromycin	E-15	15µg, 50 disks, 1 cartridge	230793
Ethambutol (Myambutol™)	EM-25	25µg, 50 disks, 1 cartridge	231575
Ethambutol (Myambutol™)	EM-50	50µg, 50 disks, 1 cartridge	231576
Fosfomycin	FOS200	200µg, 50 disks, 1 cartridge	231709
Fluconazole	FCA-25	25µg, 50 disks, 1 cartridge	232201
Furazolidone	FX-100	100µg, 50 disks, 1 cartridge	230809
Gatifloxacin	GAT-5	5µg, 50 disks, 1 cartridge	231759
Gentamicin	GM-10	10µg, 50 disks, 1 cartridge	231227
Gentamicin	GM-120	120µg, 50 disks, 1 cartridge	231693
Imipenem	IPM-10	10µg, 50 disks, 1 cartridge	231644
Isoniazid	INH-1	1µg, 50 disks, 1 cartridge	231571
Isoniazid	INH-5	5µg, 50 disks, 1 cartridge	231572
Kanamycin	K-30	30µg, 50 disks, 1 cartridge	230829
Levofloxacin	LVX-5	5µg, 50 disks, 1 cartridge	231705
Lincomycin	L-2	2µg, 50 disks, 1 cartridge	230841
Linezolid	LZD-30	30µg, 50 disks, 1 cartridge	231761
Lomefloxacin	LOM-10	10µg, 50 disks, 1 cartridge	231685
Metronidazole	MET-80	80µg, 50 disks, 1 cartridge	231604
Mezlocillin	MZ-75	75µg, 50 disks, 1 cartridge	231614
Minocycline	MI-30	30µg, 50 disks, 1 cartridge	231250
Moxalactam	MOX-30	30µg, 50 disks, 1 cartridge	231610
Nafcillin	NF-1	1µg, 50 disks, 1 cartridge	230866
Nalidixic Acid	NA-30	30µg, 50 disks, 1 cartridge	230874
Neomycin	N-30	30µg, 50 disks, 1 cartridge	230882
Netilmicin	NET-30	30µg, 50 disks, 1 cartridge	231602
Nitrofurantoin	F/M-100	100µg, 50 disks, 1 cartridge	230801
Nitrofurantoin	F/M-300	300µg, 50 disks, 1 cartridge	231149
Norfloxacin	NOR-10	10µg, 50 disks, 1 cartridge	231646
Novobiocin	NB-30	30µg, 50 disks, 1 cartridge	230890
Ofloxacin	OFX-5	5µg, 50 disks, 1 cartridge	231671

Product Name	Code	Unit of Measure	Cat. no.
Oxacillin	OX-1	1µg, 50 disks, 1 cartridge	230906
Oxolinic Acid	OA-2	2µg, 50 disks, 1 cartridge	231566
Oxytetracycline	T-30	30µg, 50 disks, 1 cartridge	230990
P-Aminosalicylic Acid	PAS-50	50µg, 50 disks, 1 cartridge	231574
Penicillin	P-2	2µg, 50 disks, 1 cartridge	230914
Penicillin	P-10	10µg, 50 disks, 1 cartridge	230918
Piperacillin	PIP-100	100µg, 50 disks, 1 cartridge	231608
Piperacillin/Tazobactam	TZP-110	110µg, 50 disks, 1 cartridge	231691
Polymyxin B	PB-300	300µg, 50 disks, 1 cartridge	230926
Quinupristin-Dalfopristin (Synercid)	SYN15	15µg, 50 disks, 1 cartridge	231720
Rifampin	RA-5	5µg, 50 disks, 1 cartridge	231541
Rifampin	RA-25	25µg, 50 disks, 1 cartridge	231578
Sparfloxacin	SPX-5	5µg, 50 disks, 1 cartridge	231707
Spectinomycin	SPT-100	100µg, 50 disks, 1 cartridge	231636
Streptomycin	S-10	10µg, 50 disks, 1 cartridge	230942
Streptomycin	S-50	50µg, 50 disks, 1 cartridge	231570
Streptomycin	S-300	300µg, 50 disks, 1 cartridge	231694
Sulbactam/Ampicillin	SAM-20	20µg, 50 disks, 1 cartridge	231659
Sulfa/Trimethoprim	SXT	23.75µg/1.25µg, 50 disks, 1 cartridge	231536
Sulfachloropyridazine	SL-1.0	1.0mg, 50 disks, 10 cartridges	231326
Sulfadiazine	SD-0.25	0.25mg, 50 disks, 10 cartridges	231331
Sulfamethizole	TH-0.25	0.25mg, 50 disks, 10 cartridges	231345
Sulfathiazole	ST-0.25	0.25mg, 50 disks, 10 cartridges	231337
Sulfisoxazole	G-0.25	0.25mg, 50 disks, 1 cartridge	230813
Sulfisoxazole	G-2.0	2.0mg, 50 disks, 1 cartridge	230821
Sulfisoxazole	G-1.0	1.0mg, 50 disks, 10 cartridges	231297
Tetracycline	Te-5	5µg, 50 disks, 1 cartridge	230994
Tetracycline	Te-30	30µg, 50 disks, 1 cartridge	230998
Ticarillin	TIC-75	75µg, 50 disks, 1 cartridge	231618
Tigecycline	TGC-15	15µg, 50 disks, 1 cartridge	232087
Timentin (Ticarillin with Clavulanic Acid)	TIM-85	85µg, 50 disks, 1 cartridge	231648
Tobramycin	NN-10	10µg, 50 disks, 1 cartridge	231568
Trimethoprim	TMP-5	5µg, 50 disks, 1 cartridge	231600
Triple Sulfa	SSS-0.25	0.25mg, 50 disks, 1 cartridge	231018
Triple Sulfa	SSS-1.0	1.0mg, 50 disks, 1 cartridge	231022
Vancomycin	Va-5	5µg, 50 disks, 1 cartridge	231030
Vancomycin	Va-30	30µg, 50 disks, 1 cartridge	231034

# Accessory Products for AST

Calipers for measuring zones, digital readout 62379531

Battery and NIST traceable certificate included

## Dispensers

Single, metal, BD 260457


6-place, self-tamping, BD 260661

8-place, self-tamping, BD 260660


12-place, self-tamping, BD 260640



12-place Dispenser,  
Cat. no. 260640



Single Disk Dispenser



Digital Calipers

## Storage Container

Holds HardyDisks™, canisters for refrigerated storage


See through latch lid, holds approximately 1751  
45 tubes, 14 1/8" L x 7 7/8" W, 4 7/8" H"



HardyDisks™ Storage  
Container, Cat. no. 1751

## DesiView™, desiccant for disk dispensers with color indicator crystals

5gm, 10 packets/box DV10



DesiView™, Cat. no. DV10

Color indicator crystals turn from blue to pink when exposed to moisture.

## McFarland Standard, 0.5 latex equivalent, turbidity standard

16x100mm tube, 8ml ML05

## Sterile Swabs

Individually wrapped Dacron® swabs, 100/pack HD25806

## Mixer II, Vortex


6.25in x 4in x 5.1in deep, each M37615

## Wickerham Card for turbidity reading

2in x 3in, each Z08

## Zone Size Interpretation Kit

Template with decals, each 260639



Wickerham Card,  
Cat. no. Z08

# Culture Media for AST

## Haemophilus Test Medium Agar

100mm plate, 10 plates/package	G33
150mm plate, 10 plates/package	H07

## Mueller Hinton Agar

100mm plate, 10 plates/package	G45
150mm plate, 10 plates/package	H11

## Mueller Hinton Agar with 5% Sheep Blood

100mm plate, 10 plates/package	A59
150mm plate, 10 plates/package	H21

## Mueller Hinton Broth

13x100mm tube, 3ml, 20 tubes/package	R38
16x100mm tube, 5ml, 20 tubes/package	K91
16x125mm tube, 10ml, 20 tubes/package	K38
500ml polycarbonate bottle, each	U31

## BHI (Brain Heart Infusion) Broth


16x100mm tube, 5ml, 20 tubes/package	K27
--------------------------------------	-----

## Saline, 0.85%

16x100mm tube, 5ml, 20 tubes/package	K59
--------------------------------------	-----



Mueller Hinton Agar, Cat. no. H11




Mueller Hinton Broth, Cat. no. U31

Mueller Hinton Agar,  
Cat. no. H11









800.266.2222  
805.346.2760 (fax)  
[www.HardyDiagnostics.com](http://www.HardyDiagnostics.com)  
[Sales@HardyDiagnostics.com](mailto:Sales@HardyDiagnostics.com)



The Quality Management System at Hardy Diagnostics manufacturing facility is certified to ISO 13485.

ISO 13485  
FM 572526

Santa Maria, CA | Centralia, WA | Salt Lake City, UT | Phoenix, AZ | Dallas, TX | Springboro, OH | Lake City, FL | Albany, NY